


V79 COLONY FORMING ASSAY

Experiment Name : ^{137}Cs toxicity (cluster) γ - irradiation

Experiment No.: B06V79.001

Date: 12/7/01

1. Seed a 225 cm² flask containing 50 ml MEMA with 3×10^6 cells 3 days before experiment.
2. Wash cells in 225 cm² flask with 20 ml PBS, trypsinize cells with 2ml Trypsin-EDTA, and resuspend in 7 ml MEMB; ^{3 min, 37°C}
10ml Pen. Str for 1 l PBS.
3. Aspirate cells and transfer to 14 ml 17x100mm Falcon tube and vigorously pass five times through 5 cc syringe with 21 gauge needle;
4. Perform cell count by transferring 100 μl in Coulter cup containing 20 ml Isotone-II (Coulter balanced electrolyte solution);
Instrument settings: Lower Threshold (LT) - 5.0; Manometer - 100 μl .
Perform cell count in triplicates. For example, 5353 5442 5291. Average value - 5362.
Cell concentration - $5362 \times 400 \times 5 = 10,7 \times 10^6$ cells/ml.

5. Dilute to ~4,000,000 cells/ml in MEMB (final volume 11 ml). For example, 4.1 ml stock + 6.9 ml MEMB gives 4×10^6 cells/ml. *Repeat Coulter counting?*
6. Transfer 1 ml of cell suspension into ten 14 ml tubes (Falcon 352059 polypropylene culture tube, 17x100 mm) labeled 1-10 both on cap and wall *fill 10 ml level*
Roll the tubes on rocker-roller for 4 h at 37°C, 5% CO₂ *10ml* **Date/Time:** 
7. After ~4 h incubation period, remove tubes, and centrifuge at 2000 rpm RT for 10 min.
8. Decant supernatant, click tubes, vortex, resuspend in 10 ml MEMA
9. Centrifuge tubes for 10 min at 2000 rpm, RT *leave ~200 μl*
10. Decant supernatant, click tubes, vortex, transfer the cell suspension to 400 μl polypropylene microcentrifuge tubes with attached caps (VWR) using 200 μl pipetter
11. Again add 200 μl MEMA, resuspend and transfer the cell suspensions in the same polypropylene microcentrifuge tubes (Total volume ~400 μl)
12. Centrifuge tubes for 5 min at 1000 rpm, RT
13. Transfer tubes at 10.5°C for 72 h. **Date/Time:** Dec 12, 2001 / 22:00

14. After 72 h, place the tubes on the perforated plate of Rainin pipet tip box containing ice (to maintain ~ 10.5°C) Dec. 10
← 2000
15. The tubes were irradiated using Mark I irradiator (¹³⁷Cs gamma-ray), one tube at a time, while placing onto a separate Rainin pipet tip box containing ice as per the Table below

Tube #	Total Dose (Rad)	Dose rate (Rad/min)	Time (min)	Attenuat.
1	0	0	0	0
2	0	0	0	0
3	100 1500	101.4	0.98	X-10 X-5
3	200	101.4	1.97 1.27	X-10 X-5
4	300	101.4	2.95 1.90	X-10 X-5
5	400	101.4	3.94 2.53	X-10 X-5
6	500	101.4	4.93 3.16	X-10 X-5
7	750	169.6	4.42 1.99 4.75	X-5 X-2
8	1000	169.6	5.89 2.66 6.33	X-5 X-2
9	1250	169.6	7.37 3.33 7.40	X-5 X-2
10	1500		8.99 9.50	X-2

1.2
2.2
3.2
4.2
5.2
6.2 6.3
7.2 7.3
8.2 8.3
9.3 9.4
10.3 10.4

3 dishes
3 ϕ 60 mm
3 * Dose rate
X-5
3 157.95 $\frac{\text{rad}}{\text{min}}$
3
6 X-2
6 375.90
6
6
6

16. After irradiation, carefully remove the supernatant from the top, resuspend pellet in 200 μ l wash MEMA and transfer the contents to ten 14 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using pasteur pipet
17. Again add 200 μ l wash MEMA in microcentrifuge tubes, resuspend and transfer the cell suspensions to 14 ml tubes
18. Centrifuge the tubes for 7 min at 2000 rpm, 4°C
19. Labeling and preparation of dilution tubes and colony dishes
 - load 57 60 mm petri dishes with 4 ml MEMA
 - load 39 T-tubes with 4.5 ml MEMA and label them [1.2, 1.3, 1.4, 1.5] [2.2, 2.3, 2.4, 2.5] X.2, X.3, X.4, X.5 etc. (1:10 dilution) previously sterilized
most diluted
less diluted
20. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
21. Centrifuge tubes for 10 min at 2000 rpm, 4°C
22. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times

6.2 \equiv tube 6 (2×10^2 /ml dilution)
6.3 \equiv tube 6 (2×10^3 /ml dilution)

Note to (23).

5379 5074 5001

$$5000(400) = 2 \times 10^6 / \text{ml}$$

manometer - 500 μl .

Q

through 5 cc syringe with 21 gauge needle

23. Determine cell concentration by transferring 100 μ l to Coulter cup (monometer - 50 μ l)
24. Vortex tube, transfer 0.5 ml into dilution tube X.5, vortex tube X.5 and transfer 0.5 ml to tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3 and vortex tube X.3 and transfer 0.5 ml to tube X.2. Keep tubes on ice. 4.0 ml
25. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.
26. Incubate petridishes for 1 week
27. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with ethanol. Stain colonies with 0.0 5% crystal violet
28. Count colonies. There must be between 25 and 250 colonies for the flask to be a valid data point.

Note: +

(24) x.5 \rightarrow 10^6 cells in 5 ml \rightarrow 2×10^5 cells/ml . OK!
x.4 \rightarrow 2×10^4 cells/ml
x.3 \rightarrow 2×10^3 cells/ml
x.2 \rightarrow 200 cells/ml

TABLE-3

Expt. # :

Date/Time :

Tube #	Coulter count for 100 ul cell suspension 500 µl manometer 20 ul Isotone II	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]
1	500 5379 5074 5001			
2	4820 4843 4898			
3	4830 4939 4866			
4	5488 5562 5577			
5	5255 5229 5452			
6	5881 5881 5850			
7	6288 6288 6277			
8	5703 5603 5520			
9	5151 5198 5125			
10	5162 5229 5172			


Dec 26, 2001

Colonies

Tube #			SF	Vol (cc)
1.2	124, 99, 109	}	1	0
2.2	91, 98, 113		105.7	
3.2	75, 92, 84	83.6	0.79	200
4.2	125, 106, 121	117.3	1.11	100
5.2	97, 96, 101	98	0.93	200
6.2	101, 103, 124	109.3	1.03	100
7.2	96, 93, 84	91	0.85	100
8.2	46		?	
8.3	210		?	1000
9.4	↑ 400		?	
10.3	128		?	1200
			?	0.12

Cluster Cg-137

SF



Small Logarithmic
Scales x 10 to the Inch

500

1000

1500

Dose (rads)

TABLE-4

Expt # :

Date :

Tube.dilution	Colony 1	Colony 2	Colony 3	Avg Colony	SF