

7

V79 COLONY FORMING ASSAY

Experiment Name : $^3\text{HTdR}$ toxicity (cluster, 50% labeling) \pm 100 μM lindane
Investigator: R Howell

Exp. #: 7
Date: 10/22/01

Preparation of cells: confluent T25 flask split 1:5 into T225 on 10/19/01. Cells trypsinized from T225 at 1:30 pm on 10/22/01. About 70% confluent.

Serum/Lot #s:

1. Set the rocker-roller at 37°C incubator with 5% CO₂, set the Coulter Counter, wash cells (from one 225 cm² flasks, seeded with 25 x 10⁶ cells the day before) with PBS-PS, trypsinize cells with 2 ml trypsin 3 min at 37°C, resuspend in 10 ml MEMB, pool, pass five times through 10 cc syringe with 21 gauge needle, perform cell count by transferring 100 μl in Coulter cup containing 20 ml Isotone II (Coulter balanced electrolyte solution).
2. Dilute to ~2,000,000 cells/ml in MEMB [Actual count : cells/ml)
3. Transfer 1 ml of cell suspension into 20 14 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall
4. Keep the tubes in the roller for 3-4 h at 37°C, 5% CO₂ **Date/Time:** 2:10 pm
5. Prepare MEMB containing radioactivity in hood
1100 μl $^3\text{HTdR}$ (Stock : 1 $\mu\text{Ci}/\mu\text{l}$ on 10/15/01) + 4.4 ml MEMB
Manufacturer: Berkin Elmer Lot #: 3106449 Calibration: 10/15/01
NET-0272
6. After 3-4 h, remove first set of 10 test tubes from roller and add MEMB with or without radioactivity according to Table below. Also add 1 ml MEMB to each of second set of tubes.

Tube #	³ HTdR μCi/ml	Cells in MEMB (ml)	MEMB (ml)	MEMB+ ³ HTdR [200 μCi/ml] (ml)
1	0	1.0	1.0	0
2	25	1.0	0.75	0.25
3	50	1.0	0.50	0.50
4	75	1.0	0.25	0.75
5	100	1.0	0	1.0
6	0	1.0	1.0	0
7	25	1.0	0.75	0.25
8	50	1.0	0.50	0.50
9	75	1.0	0.25	0.75
10	100	1.0	0	1.0

with
lindane

Need 5ml
@ 200μCi/ml
200
5.5

1000
1000

1100.0 μCi
1.1 ml ³HTdR
+ 4 ml MEMB

7. Return test tubes to roller for 12-14 h. Date/Time: 5:15 pm
8. Next day, while test tubes are in roller label 10 gamma-tubes (13 X 100 mm VWR glass test tube)
9. After ~12-14 h incubation period, remove all tubes and centrifuge at 2000 rpm at 4°C for 10 min (precooled centrifuge). Date/Time: 10/23 9:30am
10. Remove buckets from centrifuge and carefully remove 150 μl of supernatant from tubes containing radioactivity and place in pre-labeled gamma-tubes.
11. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
12. Centrifuge tubes for ⁶10 min at 2000 rpm, 4°C
13. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
14. Centrifuge tubes for ⁶10 min at 2000 rpm, 4°C
15. Decant supernatant, click tubes, vortex, resuspend in 10 ml of wash MEMA
16. Centrifuge tubes for ⁶10 min at 2000 rpm, 4°C
17. Decant supernatant, click tubes, vortex, resuspend unlabeled cells in ⁷7 ml MEMA.
 For tubes 6-10 labeled & unlabeled resuspend in 5 ml MEMA w/ 100μM lindane
 Also need lindane MEMA for wash for tomorrow
 Sonic's stock 500X for 100μM in DMSO. Prepare 55 ml MEMA w/ 100μM lindane
0.11 ml 500X lindane + 55 ml MEMA
18. Transfer unlabeled cells to the corresponding tubes containing 2,000,000 labeled cells.
19. Wash unlabeled cell tube with 5 ml MEMA and transfer to corresponding labeled cell tube.
20. Syringe the pooled cells 5 times with 5 ml syringe with 21-G needle.
21. Centrifuge tubes for 10 min at 2000 rpm, 4°C.
22. Decant supernatant completely, click tubes, vortex.

See
p. 2a

- Because mistakenly resuspended unlabeled in 7ml, had to resuspend labeled in 3ml
- Syringed labeled cells tubes 1-5 21 G 5CC 5X

Coulter count 100µl w/ 500µl manometer

1)	1543	1499	1521	608,400/ml
2)	1551	1470	1545	608,400
3)	1058	1008	989	407333
4)	856	892	867	348666
5)	872	852	881	347333

Examined cells, looked very clumpy. Decided not to make corrections to volumes based on cell counts.

Syringed all remaining tubes (labeled, unlabeled)

- Transferred 7ml of unlabeled to labeled
- Centrif. 6 min, 2000 rpm, 4°C
- Decant thoroughly, vortex, click
- Transfer 200µl to 400µl microfuge tube
- Wash 1X with 200µl MEMA or MEMA 100µl lindane (tubes 6-10)
- transfer to 400µl tube
- Spin 5 min 1000 rpm 4°C
- Transfer tubes to 10.5°C

23. Transfer the cell suspension in polypropylene microcentrifuge tubes with attached caps (Helena Plastics, 400 μ l) using 200 μ l pipette tip.
24. Again add 200 μ l ice cold MEMA, resuspend and transfer the cell suspensions in the same polypropylene microcentrifuge tubes (Total volume ~400 μ l)

25. Centrifuge tubes for 5 min at 1000 rpm, 4°C

26. Transfer tubes at 10.5°C for 72 h.

Date/Time:

10/23/01 1:00 pm

27. Transfer ²⁰30 μ l supernatant in three sets of 7 ml scintillation vials and add ⁵9 ml liquid scintillation cocktail (Ecoscint) from 150 μ l supernatant removed earlier and count them for radioactivity

Date/Time:

28. After 72 h, carefully remove the supernatant from the top, resuspend pellet in 200 μ l wash MEMA and transfer the content to ten 12 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using Pasteur pipette

Date/Time:

29. Again add 200 μ l wash MEMA in microcentrifuge tubes, resuspend and transfer the cell suspensions in 12 ml tubes

30. Centrifuge the tubes for 10 min at 2000 rpm, ^{RT}4°C (precooled centrifuge)

31. Labeling and preparation of dilution tubes and colony dishes

- load 60 mm tissue culture dishes with 4 ml MEMA

- load 40 sterile tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3, 2.4, 2.5; X.2, X.3, X.4, X.5 etc.

32. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

33. Centrifuge tubes for 10 min at 2000 rpm, ^{RT}4°C

34. ~~Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA~~

35. ~~Centrifuge tubes for 10 min at 2000 rpm, 4°C~~

36. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 5 cc syringe with 21 gauge needle

37. Determine cell concentration by transferring 100 μ l to Coulter cup

38. Vortex tube, transfer 0.5 ml into dilution tube X.5, vortex tube X.5, transfer 0.5 ml into dilution tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3, vortex tube X.3 and transfer 0.5 ml to tube X.2 and vortex. Keep tubes on ice.

39. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.

40. Transfer ¹⁰⁰200 μ l of cell suspension (in triplicate) to 20 ml scintillation vial containing 6 ml cocktail (Ecolume)

41. Incubate tissue culture dishes for 1 week

42. Count vials for radioactivity

Date/Time:

43. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol.

* May have mixed up tubes 1-5

5 tubes w/ lindane
Should wash in presence of lindane. Dilutions then remove.

Need 110 ml wash MEMA w/ lindane
0.22 ml lindane
+ 110 ml

USER: 6 ID:HOWELL H3 PRESET TIME: 1.00 TUE 30 OCT 2001 09:34
 SAMPLE REPEAT: 1 CYCLE REPEAT: 1 SCR:N R5232:N
 H 1 AQC:N GCF:N RCM:N
 CHANNEL 1-LL: 0 UL: 400 2SIGMA: 2.00 BKG SUB: 0.00 BKG 2SIG: 0.00 LSR: 0
 DATA CALC: CPM, UNKNOWN REPLICATES: 1 NORM FACTOR:Q 1.00000
 HALF LIFE(DAYS):N

SAM	POS	CH	CPM	2SIG%	TIME	EL TIME	AVG H#	ERR
1	**	1	49.00	28.57	1.00	1.42	79.0	1M
2	**	2	65.00	24.81	1.00	3.06	79.0	
3	**	3	54.00	27.22	1.00	4.63	78.0	
4	**	4	335760.00	0.89	0.15	5.35	80.0	2M
5	**	5	340640.00	0.88	0.15	6.08	80.0	
6	**	6	335030.00	1.09	0.10	6.80	80.0	3M
7	**	7	665200.00	0.78	0.10	7.55	79.0	
8	**	8	679266.62	0.63	0.15	8.31	79.0	4M
9	**	9	678900.00	0.77	0.10	9.06	80.0	
10	**	10	1015620.00	0.63	0.10	9.84	80.0	5M
11	**	11	1006173.31	0.51	0.15	10.62	79.0	
12	**	12	1013119.94	0.51	0.15	11.41	80.0	6M
13	**	13	1355813.25	0.44	0.15	12.22	80.0	
14	**	14	1372006.62	0.44	0.15	13.03	80.0	7M
15	**	15	1346480.00	0.45	0.15	13.85	80.0	
16	**	16	8.00	70.71	1.00	15.42	75.0	8M
17	**	17	4.00	100.0	1.00	16.98	74.0	
18	**	18	8.00	70.71	1.00	18.60	75.0	9M
19	**	1	336073.31	0.89	0.15	19.37	80.0	
20	**	2	340150.00	1.08	0.10	20.10	79.0	10M
21	**	3	337946.66	0.89	0.15	20.82	80.0	
22	**	4	675446.62	0.63	0.15	21.58	80.0	11M
23	**	5	682666.62	0.62	0.15	22.33	79.0	
24	**	6	670130.00	0.77	0.10	23.10	79.0	12M
25	**	7	990613.31	0.52	0.15	23.88	78.0	
26	**	8	983046.62	0.52	0.15	24.66	79.0	13M
27	**	9	1035779.94	0.51	0.15	25.45	80.0	
28	**	10	1349900.00	0.44	0.15	26.26	80.0	14M
29	**	11	1314573.25	0.45	0.15	27.07	79.0	
30	**	12	1352226.62	0.44	0.15	27.89	80.0	15M
31	**	13	11.00	60.30	1.00	29.44	93.0	
32	**	14	11.00	60.30	1.00	31.01	93.0	16M
33	**	15	4.00	100.0	1.00	32.64	92.0	
34	**	16	98446.66	1.65	0.15	33.34	93.0	17M
35	**	17	91846.66	1.70	0.15	34.04	93.0	
36	**	18	97120.00	1.66	0.15	34.74	93.0	18M
37	**	1	197240.00	1.42	0.10	35.51	92.0	
38	**	2	175070.00	1.51	0.10	36.22	94.0	19M
39	**	3	199286.66	1.16	0.15	36.94	93.0	
40	**	4	301113.31	0.94	0.15	37.66	94.0	20M
41	**	5	276480.00	1.20	0.10	38.38	92.0	
42	**	6	287820.00	1.18	0.10	39.11	93.0	21M
43	**	7	407926.66	0.91	0.15	39.84	93.0	
44	**	8	414960.00	0.98	0.10	40.57	93.0	22M
45	**	9	424200.00	0.97	0.10	41.31	94.0	
46	**	10	4.00	100.0	1.00	42.93	92.0	

SAM	POS	CH	CPM	2SIG%	TIME	EL TIME	AVG H#	ERR
47	**	-11	1	7.00	75.59	1.00	44.51	94.0
48	**	-12	1	5.00	89.44	1.00	46.14	93.0
49	**	-13	1	105500.00	1.95	0.10	46.84	92.0
50	**	-14	1	98766.66	1.64	0.15	47.55	92.0
51	**	-15	1	96853.33	1.66	0.15	48.25	94.0
52	**	-16	1	219980.00	1.35	0.10	48.97	93.0
53	**	-17	1	209893.33	1.13	0.15	49.67	92.0
54	**	-18	1	191873.33	1.18	0.15	50.39	92.0
55	**	- 1	1	308220.00	1.14	0.10	51.17	91.0
56	**	- 2	1	309030.00	1.14	0.10	51.90	93.0
57	**	- 3	1	307220.00	1.14	0.10	52.62	92.0
58	**	- 4	1	374213.31	0.84	0.15	53.35	94.0
59	**	- 5	1	352166.66	0.87	0.15	54.08	92.0
60	**	- 6	1	373870.00	1.03	0.10	54.82	93.0
61	**	- 7	1	7.00	75.59	1.00	56.43	-2.0
62	**	- 8	1	28782.86	1.99	0.35	57.40	1.0

} 6C
 } 7C
 } 8C
 } 9C
 } 10C
 Blegd.
 Std.

Exp. 7

Break up Clusters

10/26/01

Seed for Survival, Expression

Coulter Counts

Bkgd. 4 for 500ml setting

100µl cells + 20ml Iscove II + 100µl parameter

						<u>2×10^6 cells</u>
1)	1259	1259	1318	1278.7	$2.56 \times 10^6 / \text{ml}$	0.78 ml
2)	1212	1179	1241	1210.7	$2.42 \times 10^6 / \text{ml}$	0.83 ml
3)	981	970	939	963.3	$1.93 \times 10^6 / \text{ml}$	1.04 ml
4)	1022	1088	1066	1058.7	$2.12 \times 10^6 / \text{ml}$	0.94 ml
5)	1093	1038	1056	1062.3	$2.12 \times 10^6 / \text{ml}$	0.94 ml
6)	1413	1471	1471	1451.7	$2.90 \times 10^6 / \text{ml}$	0.69 ml
7)	1111	1041	1077	1076.3	$2.15 \times 10^6 / \text{ml}$	0.93 ml
8)	1057	1078	1010	1048.3	$2.00 \times 10^6 / \text{ml}$	0.95 ml
9)	1041	957	1016	1004.7	$2.00 \times 10^6 / \text{ml}$	1.0 ml
10)	891	884	937	904.0	$1.81 \times 10^6 / \text{ml}$	1.11 ml

Exp. 7

10/29/01

Replating for Expression

- Remove five T75 flasks, remove medium by aspiration.
- Wash IX 10 ml PBS, add 1 ml trypsin, 37°C, 3min
- Add 7ml MEMA, transfer to 14 ml tube
- Syringe 5X, 21G, 5cc
- Count 100µl cells, 20µl Isotone, 100µR manometer

1)	1155	1144	1191	1163.3	2.32×10^6	0.86 ml
2)	1404	1533	1438	1458.3	2.92×10^6	0.68 ml
3)	993	971	942	968.6	1.94×10^6	1.03
4)	917	895	899	903.6	1.81×10^6	1.10
5)	746	777	764	762.3	1.52×10^6	1.31
6)	1180	1121	1285	1195.3	2.39×10^6	0.84 ml
7)	1026	932	997	985	1.97×10^6	1.02 ml
8)	939	925	905	923	1.85×10^6	1.08
9)	760	732	720	737.3	1.47×10^6	1.36
10)	698	702	706	702	1.40×10^6	1.42

Exp. 7

²
11/2/01

Replating for Expression

- Aspirate media from T150, wash 1X PBS, 2 ml trypsin, 3min 37°C, resuspend in 7ml MEMA, transfer to 14 ml tube
- Syringe 5X, 21G, 5cc
- Coulter Count 500mA 5.0 Lower Thresh
100µl cells, 20µl Isotone II, 100µl macrometer

Bkg. 20

1)	2383	2564	2471	2472.7	$4.95 \times 10^6 / \text{ml}$	0.40 ml
2)	2899	2843	2909	2883.7	$5.77 \times 10^6 / \text{ml}$	0.35
3)	1808	1840	1922	1856.7	$3.71 \times 10^6 / \text{ml}$	0.54
4)	1761	1852	1789	1800.7	$3.60 \times 10^6 / \text{ml}$	0.55
5)	1777	1721	1768	1755.3	$3.51 \times 10^6 / \text{ml}$	0.57 ml
6)	3059	3060	3081	3100	$6.20 \times 10^6 / \text{ml}$	0.32
7)	2791	2736	2752	2759.7	$5.52 \times 10^6 / \text{ml}$	0.36
8)	2193	2202	2202	2199	$4.39 \times 10^6 / \text{ml}$	0.45
9)	1754	1799	1716	1756	$3.51 \times 10^6 / \text{ml}$	0.57
10)	1489	1471	1509	1489.7	$2.98 \times 10^6 / \text{ml}$	0.67

- Transfer appropriate volume of cells to T150 containing 25ml MEMA

Stain colonies with 0.05% crystal violet

44. Count colonies. There must be between 25 and 250 colonies for the dish to be a valid data point.

Tube	Dilution	Colony Counts

VT9 Exp. 7

11/5/01

Processing for Mutant Selection

- Load 100 P60 w/ MEMA 10ml * Careful
- Load 30 P60 w/ 4ml MEMA
- Load 40 dilution tubes with 4.5ml wash MEMA

- Aspirate medium from T75
- Wash 1X 10ml PBS
- Trypsinize 2ml, 3 min, 37°C, hit flasks
- Resusp in 8ml MEMA (wash) & transfer to 14ml tube
- Spin 5 min 2000 rpm, RT
- Decant, resuspend in 5ml wash MEMA
- Syringe 5X 21G 5cc

Counter Counts (100µl cells, 20µl Isotone, LT 5:10, 100µl manometer)

to get vol. for 2x10⁶ cells

1)	3438	3571	3566	3525	7.05 x 10 ⁶ /ml	0.142 cells + 0.358 MEMA
2)	3715	3784	3761	3753.3	7.50 x 10 ⁶ /ml	0.133 cells + 0.367 MEMA
3)						
4)						
5)						
6)						
7)						
8)						
9)						
10)						

V79 Exp 7.

Coulter Counts

					Cells/ml	5ml @ 2×10^6 /ml for	
						ml cells	ml MEM
1)	3438	3571	3566	3525	7.05×10^6 /ml	1.41	3.59
2)	3715	3784	3761	3753.3	7.51×10^6 /ml	1.33	3.67
3)	2164	1991	1933				
4)	1749	1714					
				Didn't spin down & resusp in 5ml			
3)	3733	3756	3801	3763.3	7.52×10^6	1.33	3.67
4)	3391	3243	3336	3323.3	6.64×10^6	1.50	3.50
5)	3322	3294	3405	3340.3	6.68×10^6	1.50	3.50
6)	3770	3693	3863	3775.3	7.55×10^6	1.32	3.68
7)	3675	3410	3525	3536.7	7.07×10^6	1.41	3.59
8)	3440	3353	3605	3466	6.93×10^6	1.44	3.56
9)	3503	3468	3397	3456	6.91×10^6	1.45	3.55
10)	3538	3371	3465	3458	6.92×10^6	1.45	3.55

- Transfer designated ml cells + ml MEMA to fresh 14 ml tube to get 2×10^6 cells/ml
- Serially dilute 0.5 ml of above to 4.5 ml wash MEMA & repeat 3X
- Seed P60's with last dilution (1ml)
- Seed P100's with 100ml of 2×10^6 cells/ml dilution

- Normally add GTG in 3 ml - need to concentrate it 10X.

- Need 30ml of GTG @ $\frac{500 \text{ Mg/ml}}{3} = 5000 \text{ Mg}$

$$\text{Stock} = \frac{5 \text{ mg}}{2 \text{ ml}} = 2.5 \frac{\text{mg}}{\text{ml}} = 2500 \text{ Mg/ml}$$

Need ~~2ml~~ 2 ml stock + 28 ml MEMA yields

$$30 \text{ ml @ } 167 \frac{\text{Mg}}{\text{ml}}$$

So add 1ml to each of 3 ~~10~~ 2.5 mg sigma vials

- Add 0.3 ml of GTG - medium to each P100