

V79 COLONY FORMING ASSAY

Experiment Name : $^{125}\text{IUdR}$ + 50-200 ug/ml MEA; Exp. # : 2; Investigator: A.Bishayee
 Date: ~~05/14/98~~ 06/15/98.

1. Set the rocker-roller at 37°C incubator, set the Coulter Counter, wash cells (from 75 cm² flask, subcultured 1:2, 24h before) with PBS, trypsinize cells, resuspend in 7 ml MEMB, pass five times through 3 cc syringe with 21 gauge needle, perform cell count by transferring 100 ul in Coulter cup containing 20 ml isotone (Coulter balanced electrolyte solution)
2. Dilute to ~4,00,000 cells/ml in MEMB (final volume 11 ml) [Actual count : 426800 cells/ml]
3. Transfer 1 ml of cell suspension into ten 12 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall
4. Roll the tubes for 3-4 h at 37°C, 5% CO₂ Date/Time: 4-15 p.m.
5. Prepare MEMB containing radioactivity in hood
3.4 µl $^{125}\text{IUdR}$ (prepared on 04/09/98) + 3 ml MEMB
6. After 3-4 h, remove test tubes from roller and add MEMB with or without radioactivity according to Table below. Date/Time: 8/15/98; 8-00 p.m.

| Tube # | $^{125}\text{IUdR}$ uCi/ml | Cells in MEMB (ml) | MEMB (ml) | MEMB+ $^{125}\text{IUdR}$ (ml) [1.0 uCi/ml] | MEMA+ MEA (200 ug/ ml) (ml) | MEMA (ml) | MEA Conc. (ug/ml) |
|--------|-------------------------------|--------------------------|--------------|---|---|--------------|-------------------------|
| 1 | 0 | 1.0 | 1.0 | 0 | 0 | 2.0 | 0 |
| 2 | 0 | 1.0 | 1.0 | 0 | 0.5 | 1.5 | 50 |
| 3 | 0 | 1.0 | 1.0 | 0 | 1 | 1 | 100 |
| 4 | 0 | 1.0 | 1.0 | 0 | 1.5 | 0.5 | 150 |
| 5 | 0 | 1.0 | 1.0 | 0 | 2.0 | 0 | 200 |
| 6 | 0.02 | 1.0 | 0.96 | 0.04 | 0 | 2.0 | 0 |
| 7 | 0.02 | 1.0 | 0.96 | 0.04 | 0.5 | 1.5 | 50 |
| 8 | 0.02 | 1.0 | 0.96 | 0.04 | 1 | 1 | 100 |
| 9 | 0.02 | 1.0 | 0.96 | 0.04 | 1.5 | 0.5 | 150 |
| 10 | 0.02 | 1.0 | 0.96 | 0.04 | 2.0 | 0 | 200 |

7. Return test tubes to roller for 12 h, increase the elevation angle of the roller.

Date/Time: 06/15/98; 8-15 p.m.

8. While test tubes are rolling label 40 (4x10) gamma-tubes (13 X 100 mm VWR glass test tube)
9. After ~12 h incubation period, remove tubes and centrifuge at 2000 rpm at 4°C for 10 min
(precooled centrifuge). Date/Time: 06/16/98; 9-10 a.m.
10. During centrifugation, move roller to 10°C and obtain ice $MEMA = 10.89 \text{ ml}$
11. Prepare 11 ml of 200 ug/ml MEA in MEMA, put on ice $MEA = 0.11 \text{ ml}$
(20 mg/ml)
12. Remove buckets from centrifuge and carefully remove 100 µl of supernatant and place in pre-labeled gamma-tube.
13. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
14. Centrifuge tubes for 10 min at 2000 rpm, 4°C
15. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
16. Centrifuge tubes for 10 min at 2000 rpm, 4°C
17. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
18. Centrifuge tubes for 10 min at 2000 rpm, 4°C
19. Decant supernatant, click tubes, vortex, resuspend in 2 ml ice cold MEMA containing 0 and/or 200 ug/ml MEA as per Table. Keep on ice!
20. Transfer tubes to roller at 10°C for 72 h. Date/Time: 06/16/98; 11-00 a.m.
21. Transfer 10 ul supernatant in three sets of tubes containing small pieces of tissue paper from 100 ul supernatant removed earlier and count them for radioactivity
Date/Time: 06/16/98; 12-00 noon
21. After 72 h, remove tubes and place on ice, add 8 ml ice cold wash MEMA.
Date/Time: 06/19/98; 9-20 a.m.
22. Centrifuge tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)
23. Labeling and preparation of dilution tubes and colony dishes
 - load 57 60 mm petri dishes with 4 ml MEMA
 - load 30 T-tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 2.2, 2.3, 2.4, X.2, X.3, X.4, etc.
24. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
25. Centrifuge tubes for 10 min at 2000 rpm, 4°C
26. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
27. Centrifuge tubes for 10 min at 2000 rpm, 4°C
28. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 3 cc syringe with 21 gauge needle
29. Determine cell concentration by transferring 100 µl to Coulter cup
30. Vortex tube, transfer 0.5 ml into dilution tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3, vortex tube X.3 and transfer 0.5 ml to tube X.2 and vortex. Keep tubes on ice.
31. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.

32. Transfer 300 μ l of cell suspension (in triplicate) to gamma tubes for each tube
33. Incubate petridishes for 1 week
34. Count gamma tubes for radioactivity Date/Time : 06/19/98; 4-30 p.m .
35. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol.
 Stain colonies with crystal violet
36. Count colonies (50 or more cells). There must be between 25 and 250 colonies for the flask to be a valid data point.

Expt # 2

06/15/98

$$\begin{aligned}
 \text{Initial Cell count} &= 8693, 8716, 8638 \\
 \text{Avg. Cell count} &= 8682.3 \\
 \text{Cell conc.} &= 8682.3 \times 400 \\
 &= 3472933 \text{ cells/ml}
 \end{aligned}$$

For dilution,

$$\begin{aligned}
 \text{Vol of cell suspension taken} &= \frac{4000000}{3472933} \\
 &= 1.15 \text{ ml}
 \end{aligned}$$

Take 1.3 ml of Cello + 9.7 ml MEMB = 11 ml.

After dilution,

$$\begin{aligned}
 \text{Final Cell count} &= 1121, 1070, 1010 \\
 \text{Avg. Cell count} &= 1067 \\
 \text{Cell conc.} &= 1067 \times 400 \\
 &= 426800 \text{ Cello/ml.}
 \end{aligned}$$

Expt #2

06/15/98

Preparation of ^{125}I Udr in MEMB

Prepare 3 ml of 1 $\mu\text{Ci}/\text{ml}$ = 3 μCi required

Stock

$$\text{on } 4/9/98 = 1.93 \mu\text{Ci}/\text{ml}$$

$$\begin{aligned} \text{on } 6/15/98 &= 1.93 \times 0.4588 \\ &= 0.88 \mu\text{Ci}/\text{ml}. \end{aligned}$$

$$\text{Stock required} = \frac{3}{0.88} = 3.4 \mu\text{l}.$$

① Take 3 ml of MEMB

② Take 3.4 μl Stock ^{125}I Udr

TABLE-1

Expt. #: 2

Date/Time: 06/16/98; 12-00 noon

| Tube # | Medium count for 10 ul (cpm) | Avg. cpm | dpm [cpm/0.7056] | $\mu\text{Ci/ml (A}_1)$ on counting [dpm/22200] | $\mu\text{Ci/ml (A}_2)$ on addition [A ₁ e ^{-λt}] | |
|--------|---------------------------------|----------|---------------------|---|--|-------|
| 1 | 0, 1, 1 | 0 | 0 | 0 | 0 | |
| 2 | 0, 0, 0 | 0 | 0 | 0 | 0 | |
| 3 | 0, 1, 0 | 0 | 0 | 0 | 0 | |
| 4 | 0, 0, 1 | 0 | 0 | 0 | 0 | |
| 5 | 1, 1, 0 | 0 | 0 | 0 | 0 | |
| 6 | 270, 272, 260, | 267.3 | 378.8 | 0.0017 | 0.0017 | 0.017 |
| 7 | 246, 216, 280 | 247.3 | 350.5 | 0.0015 | 0.0015 | 0.015 |
| 8 | 293, 255, 320 | 289.3 | 410.0 | 0.0018 | 0.0018 | 0.018 |
| 9 | 252, 276, 244 | 257.3 | 364.7 | 0.016 | 0.016 | 0.016 |
| 10 | 260, 258, 236 | 251.3 | 356.1 | 0.016 | 0.016 | 0.016 |

06/15/98; 8:05 P.M.

12h + 4h

= 16h

$$e^{-\lambda t} = e^{-\frac{0.693 \times 16}{1440}}$$

$$= 0.9923$$

TABLE-2

Expt. # :

Date/Time : 06/19/98; 4-30 p.m.

| Tube # | Radioactivity for 300 ul cell suspension (cpm) | Avg. cpm | dpm [cpm/0.7056] | $\mu\text{Ci/ml (A)}$ on counting [dpm/666000] | $\mu\text{Ci/ml (A}_0)$ after 12 h incubation [$A/e^{-\lambda t}$] |
|--------|--|----------|------------------|--|--|
| 1 | 0, 1, 1 | 0 | 0 | 0 | 0 |
| 2 | 0, 1, 0 | 0 | 0 | 0 | 0 |
| 3 | 2, 1, 0 | 0 | 0 | 0 | 0 |
| 4 | 2, 2, 0 | 0 | 0 | 0 | 0 |
| 5 | 0, 1, 2 | 0 | 0 | 0 | 0 |
| 6 | 1938, 1977, 2002 | 1972.3 | 2795.2 | 0.00419 | 0.00436 |
| 7 | 2145, 2140, 2030 | 2095.6 | 2970.0 | 0.00445 | 0.00463 |
| 8 | 2183, 2114, 2002 | 2099.6 | 2975.7 | 0.00446 | 0.00464 |
| 9 | 1998, 1942, 1995 | 1978.3 | 2803.7 | 0.00420 | 0.00437 |
| 10 | 1854, 1967, 1822 | 1881 | 2665.8 | 0.00400 | 0.00415 |

06/16/98; 9-00 a.m.

$$= e^{-\lambda t} = e^{-\frac{0.693 \times 79.5}{1440}}$$

$$= 0.9624$$

$$72h + 7.5 = 79.5$$

TABLE-3

Expt. # : 2

Date/Time : 06/19/98; 10-30 a.m

| Tube # | Coulter count for 100 ul cell suspension | Avg. count | Cells/ml [Avg. count x 400] | pCi/cell [uCi/ml x 10 ⁶ Cells/ml] |
|--------|--|------------|-----------------------------------|--|
| 1 | 753, 759, 734 | 748.6 | 299466 | 0 |
| 2 | 769, 780, 752 | 767 | 306800 | 0 |
| 3 | 832, 856, 842 | 843.3 | 337333 | 0 |
| 4 | 861, 868, 867 | 865.3 | 346133 | 0 |
| 5 | 766, 708, 704 | 726 | 290400 | 0 |
| 6 | 872, 862, 869 | 867.6 | 347066 | 0.0125 |
| 7 | 869, 875, 892 | 878.6 | 351466 | 0.0131 |
| 8 | 855, 862, 872 | 863 | 345200 | 0.0134 |
| 9 | 924, 908, 904 | 912 | 364800 | 0.0119 |
| 10 | 769, 777, 748 | 764.6 | 305866 | 0.0135 |

TABLE-4

Expt. #: 2

Date: 06/26/98

Colony Counts and Survival Fraction

| Tube.dilution | Colony 1 | Colony 2 | Colony 3 | Avg Colony for x.2 | SF |
|---------------|----------|----------|----------|-----------------------|--------|
| 1.2 | 110 | 119 | 105 | 113.3 | — |
| 2.2 | 108 | 99 | 105 | 104 | 0.9179 |
| 3.2 | 90 | 97 | 106 | 97.6 | 0.8620 |
| 4.2 | 80 | 66 | 72 | 72.6 | 0.6413 |
| 5.2 | 70 | 60 | 64 | 64.6 | 0.5707 |
| 6.3 | 59 | 62 | 57 | 59.3 | 0.0523 |
| 7.3 | 70 | 75 | 65 | 7.0 | 0.0673 |
| 8.3 | 225 | 220 | 230 | 22.5 | 0.2305 |
| 9.3 | 204 | 209 | 199 | 20.4 | 0.2809 |
| 10.3 | 71 | 76 | 66 | 7.1 | 0.1099 |

Conc. of MEA
($\mu\text{g/ml}$)

DMF

50

1.14

100

2.12

150

2.23

200

1.42

Expt # 2

