

V79 COLONY FORMING ASSAY

Experiment Name : $^{125}\text{IUdR}$ + MEA; **Exp. # :** 2; **Investigator:** A. Bishayee

Date: 04/27/98

1. Set the rocker-roller at 37°C incubator, set the Coulter Counter, wash cells (from 75 cm² flusk, subcultured 1:2, 24h before) with PBS, trypsinize cells, resuspend in 7 ml MEMB, pass five times through 3 cc syringe with 21 gauge needle, perform cell count by transferring 100 ul in Coulter cup containing 20 ml isotone (Coulter balanced electrolyte solution)
2. Dilute to ~4,00,000 cells/ml in MEMB (final volume 11 ml) [Actual count : 4, 22, 133 cells/ml]
3. Transfer 1 ml of cell suspension into ten 12 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall
4. Roll the tubes for 3-4 h at 37°C, 5% CO₂ **Date/Time:** 04/27/98; 3-45 p.m.
5. Prepare MEMB containing radioactivity in hood
3.2 μl $^{125}\text{IUdR}$ (prepared on 04/09/98) + 5 ml MEMB
6. After 3-4 h, remove test tubes from roller and add MEMB with or without radioactivity according to Table below. **Date/Time:** 04/27/98; 7-15 p.m.

Tube #	$^{125}\text{IUdR}$ uCi/ml	Cells in MEMB (ml)	MEMB (ml)	MEMB+ $^{125}\text{IUdR}$ (ml) [1.0 uCi/ml]	MEMA+ MEA 100 ug/ml (1.3 mM)	MEMA	
1	0	1.0	1.0	0	2.0	0	
2	0	1.0	1.0	0	2.0	0	
3	0.01	1.0	0.98	0.02	2.0	0	
4	0.02	1.0	0.96	0.04	2.0	0	
5	0.03	1.0	0.94	0.06	2.0	0	
6	0	1.0	1.0	0	0	2.0	
7	0	1.0	1.0	0	0	2.0	
8	0.01	1.0	0.98	0.02	0	2.0	
9	0.02	1.0	0.96	0.04	0	2.0	
10	0.03	1.0	0.94	0.06	0	2.0	

7. Return test tubes to roller for 12 h, increase the elevation angle of the roller.

Date/Time: 04/27/98; 7-30 p.m.

1.2

2.2

3.2, 3.3

4.2, 4.3

5.3, 5.4

6.2

7.2

8.2, 8.3

9.2, 9.3, 9.4

10.3, 10.4

8. While test tubes are rolling label 40 (4x10) gamma-tubes (13 X 100 mm VWR glass test tube)

9. After ~12 h incubation period, remove tubes and centrifuge at 2000 rpm at 4°C for 10 min

(precooled centrifuge).

Date/Time: 04/28/98 (MEA (20 ng/ml) = 55 μl

10. During centrifugation, move roller to 10°C and obtain ice ^{9-10 a.m.}

MEMA = 10.945 ml

11. Prepare 11 ml of sterile MEA (100 ug/ml) in MEMA, put on ice

= 10.945 ml

12. Remove buckets from centrifuge and carefully remove 100 μl of supernatant and place in pre-labeled gamma-tube.

13. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

14. Centrifuge tubes for 10 min at 2000 rpm, 4°C

15. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

16. Centrifuge tubes for 10 min at 2000 rpm, 4°C

17. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

18. Centrifuge tubes for 10 min at 2000 rpm, 4°C

19. Decant supernatant, click tubes, vortex, resuspend in 2 ml ice cold MEMA containing 0 or 100 ug/ml of MEA as per Table. Keep on ice!

20. Transfer tubes to roller at 10°C for 72 h.

Date/Time: 04/28/98; 11-00 a.m.

21. Transfer 10 ul supernatant in three sets of tubes containing small pieces of tissue paper from 100 ul supernatant removed earlier and count them for radioactivity

Date/Time: 04/28/98;

21. After 72 h, remove tubes and place on ice, add 8 ml ice cold wash MEMA.

Date/Time: 05/01/98; 9-10 a.m.

22. Centrifuge tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)

23. Labeling and preparation of dilution tubes and colony dishes

- load 57 60 mm petri dishes with 4 ml MEMA

- load 30 T-tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 2.2, 2.3, 2.4, X.2, X.3, X.4, etc.

24. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

25. Centrifuge tubes for 10 min at 2000 rpm, 4°C

26. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

27. Centrifuge tubes for 10 min at 2000 rpm, 4°C

28. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 3 cc syringe with 21 gauge needle

29. Determine cell concentration by transferring 100 μl to Coulter cup

30. Vortex tube, transfer 0.5 ml into dilution tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3, vortex tube X.3 and transfer 0.5 ml to tube X.2 and vortex. Keep tubes on ice.

31. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.

32. Transfer 300 μ l of cell suspension (in triplicate) to gamma tubes for each tube
33. Incubate petridishes for 1 week
34. Count gamma tubes for radioactivity **Date/Time : 05/04/98; 12-00 noon**
35. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol.
Stain colonies with crystal violet
36. Count colonies (50 or more cells). There must be between 25 and 250 colonies for the dish to be a valid data point.

Expt # 2

04/27/98

Initial cell count = 6495, 6741, 6880, 6741
Avg. cell count = 6787.3
Cell conc. = 400×6787.3
= 2714933.3 cells/ml

For dilution

$$\text{Vol. of cell suspension required} = \frac{4400000}{2714933.3} = 1.62$$

Take 1.62 ml of cells + 9.38 ml MEMB = 11 ml.

After dilution,

Final count = 1127, 1085, 954
Avg. cell count = 1055.3
Cell conc = 400×1055.3
= 4,22,133 cells/ml

Expt #2

04/27/98

Preparation of ^{125}I IWR in MEMB.

Prepare 5 ml of $1 \mu\text{Ci}/\text{ml}$ ^{125}I IWR = 5 μCi required.

Stock on 04/09/98 $1.93 \mu\text{Ci}/\text{ml}$

on 04/27/98 $1.93 \times 0.812 \mu\text{Ci}/\text{ml}$
 $= 1.57 \mu\text{Ci}/\text{ml}$.

Stock required = $\frac{5}{1.57} = 3.19 \text{ ml}$.

- ① Take 5 ml of MEMB
- ② Add 3.2 ml of stock ^{125}I IWR

Expt Name : ¹²⁵IUR + 1.3 mM MEA

Expt # : 2

Date / Time : 04/28/98 ; 12-00

10 µl onto time

1	1.00	Blank	20	17.5	234
2	1.00		15		229
3	1.00	1M	28		249
4	1.00		17		252
5	1.00	2M	23		239
6	1.00		25		247
7	1.00	3M	18		232
8	1.00		25		247
9	1.00	4M	172		431
10	1.00		195		412
11	1.00	5M	197		384
12	1.00		342		571
13	1.00	6M	317		512
14	1.00		314		526
15	1.00	7M	451		670
16	1.00		487		727
17	1.00	8M	455		656
18	1.00		22		201
19	1.00	9M	17		219
20	1.00		26		221
21	1.00	10M	17		226
22	1.00		24		228
23	1.00	Blank	21		255
24	1.00		190		435
25	1.00	Blank	182		390
26	1.00		189		382
27	1.00	Blank	335		566
28	1.00		309		511
29	1.00	Blank	348		578
30	1.00		479		687
31	1.00	Blank	511		739
32	1.00		502		706
33	1.27	Blank	8		68

TABLE-1

Expt. # : 2

Date/Time : 04/28/98; 12-00

Tube #	Medium count for 10 ul (cpm)	Avg. cpm	dpm ^{0.7056} [cpm/0.142]	μCi/ml (A) on counting [dpm/22200]	μCi/ml (A ₀) on addition [A _t e ^{-λt}]
1	11, 0, 6	0	0	0	0
2	8, 1, 8	0	0	0	0
3	155, 178, 180	171	242.3	0.0109	0.0110
4	325, 300, 297	307.3	435.5	0.0196	0.0197
5	434, 470, 438	447.3	633.9	0.0285	0.0287
6	5, 0, 9	0	0	0	0
7	0, 7, 4	0	0	0	0
8	173, 165, 172	170	240.9	0.0108	0.0109
9	318, 292, 331	313.6	444.5	0.0200	0.0201
10	462, 494, 485	480.3	680.7	0.0306	0.0309

04/27/98; 7-30 p.m.

$$= e^{-\lambda t}$$

$$= e^{-\frac{0.693 \times 16.5}{1440}}$$

$$= e^{-0.9920}$$

$$12h + 4.5h$$

$$= 16.5h$$

Expt Name : ¹²⁵IodR+ 1.3mM MEA

Expt # : 2

Date : 05/04/98

Time : 12-00 noon.

POWER-UP. MODEL 9455 VERS. 1.0 PRINTER UPDATE.

300ul cells.

POWER-UP. MODEL 9455 VERS. 1.0 PRINTER UPDATE.

1	1.00		{ 31	Back.	248
2	1.00	26	{ 22		263
3	1.00		{ 18		196
4	1.00		{ 13	1c	238
5	1.00		{ 24		230
6	1.00		{ 25		223
7	1.00		{ 20	2c	228
8	1.00		{ 30		221
9	1.00		{ 558		774
10	1.00		{ 515	3c	745
11	1.00		{ 538		753
12	1.00		{ 910		1115
13	1.00		{ 878	4c	1088
14	1.00		{ 894		1129
15	1.00		{ 1635		1857
16	1.00		{ 1628	5c	1854
17	1.00		{ 1520		1746
18	1.00		{ 21	6c	238
19	1.00		{ 25		229
20	1.00		{ 25		228
21	1.00		{ 18	7c	251
22	1.00		{ 18		210
23	1.00		{ 27		222
24	1.00		{ 569		754
25	1.00		{ 595	8c	796
26	1.00		{ 642		853
27	1.00		{ 871	9c	1089
28	1.00		{ 904		1098
29	1.00		{ 880		1114
30	1.00		{ 1423		1614
31	1.00		{ 1480	10c	1695
32	1.00		{ 1564		1779

TABLE-2

Expt. # : 2

Date/Time : 05/04/98; 12-00 noon

Tube #	Radioactivity for 300 ul cell suspension (cpm)	Avg. cpm	dpm [cpm/0.142] 0.7056	μ Ci/ml (A) on counting [dpm/666000]	μ Ci/ml (A ₀) after 12 h incubation [A _t e ^{-λt}]
1	-8, -13, -2	0	0	0	0
2	-1, -6, 4				
3	532, 489, 512	511	724.20	0.00108	0.00116
4	884, 852, 868	868	1230.15	0.00184	0.00198
5	1609, 1602, 1494	1568.3	2222.69	0.00333	0.00358
6	-5, -1, -1	0	0	0	0
7	-8, -8, 1	0	0	0	0
8	543, 569, 616	576	816.32	0.00122	0.00131
9	845, 878, 854	859	1217.40	0.00182	0.00196
10	1397, 1454, 1538	1463	2073.4	0.00311	0.00334

$$e^{-\lambda t}$$

$$= e^{-\frac{0.693 \times 147}{1440}}$$

$$= 0.9317$$

04/20/98; 9-10 a.m.

$$6 \times 24 \text{ h} + 3$$

$$= 147 \text{ h.}$$

TABLE-3

Expt. # : 2

Date/Time : 05/01/98; 11-15 a.m.

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]
1	589, 529, 512	543.3	217333.3	0
2	564, 560, 585	569.6	227866.6	0
3	691, 627, 605, 593	608.3	243333.3	0.00476
4	518, 461, 460	479.6	191866.6	0.01031
5	695, 667, 695	685.6	274266.6	0.01305 0.06756
6	730, 721, 749	733.3	293333.3	0
7.	621, 599, 613	611	244400	0
8	672, 687, 666	675	270000	0.00485
9	628, 576, 523, 533	544	217600	0.009007
10	458, 421, 475	451.3	180533.3	0.0185007

TABLE-4

Expt #: 2

Date: 05/08/98

Colony Counts and Survival Fraction

Tube.dilution	Colony 1	Colony 2	Colony 3	Avg Colony	SF
1:2	58	54	56	} 60.5	
2:2	69	61	65		
3:2	39	33	37	36.33	0.6005
4:3	82	97	91	9.0	0.1487
5:3	58	45	46	4.96	0.0820
6:2	70	75	78	} 71	
7:2	69	68	66		
8:2	23	29	22	24.66	0.3474
9:3	38	40	30	3.6	0.0507
10:4	27	25	36	0.2933	0.0041

Expt # 2

○ — ○ with MEA
⊗ — ⊗ without MEA

DMF = 1.8

