

15

1

## V79 COLONY FORMING ASSAY

**Experiment:** 100 cluster conformation, Cs-137 rays);

**Exp.**

**Investigator:** M. Lenarczyk

**Date:** April 2, 2001

### ***Procedure to get V79 cells with low background for HPRT- mutants.***

Cells were propagated in normal culture medium upto 70-80 confluency.

Then they were split into 2 T175 flask and they were under HAT (1x) condition (from HAT 100x stock solution)

Next day they the cells were trypsinized @ split in same 2 T175 flasks and were cultured under HAT condition

After two days under HAT culture condition cells were split into 2 T175 flasks and were cultured for next two days in normal medium without HAT.

Day before experiment cells were re-plated (app.  $10 \times 10^6$ ) into each of two T175 ( from app. 80% confluency =  $30.5 \times 10^6$ ) (flask A) and  $31.5 \times 10^6$ ) (flask C).

Next day cells were washed with HBSS, trypsinized and counted

**Date:** April 02, 2001

Flask A = count -  $162 \times 10^4$  cells x 20 ml =  $32.4 \times 10^6$  cells (app. 80-90% confluency)

Flask C = count -  $168 \times 10^4$  cells x 20 ml =  $33.6 \times 10^6$  cells (app 80-90 % confluency)

### ***Experimental procedure***

1. Set the rocker-roller at 37°C incubator with 5% CO<sub>2</sub>, pool, pass 5x through 3 cc syringe with 21 gauge needle,
2. Dilute to ~4,000,000 cells/ml in MEMB [Actual count : 4 000 000 cells/ml]
  - a) cells were counted using hemacytometer
  - b) flask A + flask C =  $(32.4 + 33.6) \times 10^6$  cells =  $66 \times 10^6$  cells
    1. re-suspend  $66 \times 10^6$  cells in 10 ml S-MEM
    2. take 9.09 ml (app.  $60 \times 10^6$  cells) and add 5.91 ml S-MEM =  $60 \times 10^6$  / 15 ml =  $4 \times 10^6$ /ml
    3. count the cells : counts is  $(306 + 322) / 2 \times 10^4$  cells =  $3.14 \times 10^6$
    4. take 1.27 ml ( $4 \times 10^6$  cells) and transfer them into Falcon polypropylene test tube, (7x100 mm) labeled 1-10 both on cap and wall 12 ml Falcon tube
    5. add 0.73 ml S-MEM into each tube to make final volume = 2.00 ml

Sample	Cs-137 (Gy)	Temp. (10.5° C)
1	control	100% cluster
2	control	100% cluster
3	Gy	100% cluster
4	Gy	100% cluster
5	Gy	100% cluster
6	Gy	100% cluster
7	Gy	100% cluster
8	Gy	100% cluster
9	Gy	100% cluster
10	Gy	100% cluster
11	Gy	control after ROLL

\* Before I put the cells in clusters I also plated them for survival to check "background" PE

$74 \times 10^4$  cells  $\frac{\text{cluster } 10^2}{5}$   $74 \times 10^3$   $\frac{\text{cluster } 10^2}{5}$   $74 \times 10^2 / \text{ml}$

0.203 ml of  $74 \times 10^2 / \text{ml} = 1500$  cells

0.203 ml + 29.8 ml MEMA = 1500 cells / 30 ml

$\swarrow$  2 ml / P60  
 $\searrow$  1 ml / P60  
 ↓  
 100 cells / P60 + 2 ml MEMA × 3  
 200 cells / P60 × 3  
 ↓  
 37°C, 7 days

\* \* \* →

3. Keep the tubes in the roller overnight at 37°C, 5% CO<sub>2</sub> Date/Time: April 02, 2001 / 18:00
4. After overnight incubation period (      hrs) remove tubes from incubator, and centrifuge at 2000 rpm at 4°C for 10 min (precooled centrifuge).  
Date/Time: April 3, 2001 : 10:30
5. Decant supernatant, click tubes, vortex, re-suspend in 10 ml wash MEMA
6. Centrifuge tubes for 10 min at 2000 rpm, 4°C
7. Decant supernatant, click tubes, vortex, re-suspend in 10 ml wash MEMA
8. Centrifuge tubes for 10 min at 2000 rpm, 4°C
9. Decant supernatant, click tubes, vortex, re-suspend in 10 ml wash MEMA
10. Centrifuge tubes for 10 min at 2000 rpm, 4°C
11. Decant supernatant, click tubes, vortex, re-suspend in 7 ml of MEMA (**culture medium**)
12. Centrifuge tubes for 10 min at 2000 rpm, 4°C
13. Decant supernatant, click tubes, vortex, transfer the cell suspension in polypropylene micro-centrifuge tubes Helena Plastics, (400 µl) using 200 µl pipet tips
14. Again add 200 µl ice cold MEMA, re-suspend and transfer the cell suspensions in the same polypropylene micro-centrifuge tubes (Total volume ~400 ul)
15. Centrifuge tubes for 5 min at 1000 rpm, 4°C
16. Transfer tubes at 10.5°C for 72 h. Date/Time: April 3, 2001 /

17. After 72 hrs remove cluster from 10.5 C and transfer them in plastic container contains ice.  
The cluster must be little bit above the ice.
18. Irradiate cluster with gamma rays, (Cs-137, acute exposure). Be sure that only one cluster is irradiated at time.

**Date/Time:**

19. After irradiation, carefully remove the supernatant from the top, re-re-suspend pellet in 200  $\mu$ l wash MEMA and transfer the content to ten 12 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using pasteur pipet
20. Again add 200  $\mu$ l wash MEMA in micro-centrifuge tubes, re-suspend and transfer the cell suspensions in 12 ml tubes
21. Centrifuge the tubes for 10 min at 2000 rpm, 4°C (*precooled centrifuge*)
22. Labeling and preparation of dilution tubes and colony dishes
  - load 66, 60 mm petri dishes with 4 ml MEMA
  - load 40 sterile tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3, 2.4, 2.5; X.2, X.3, X.4, X.5 etc.
23. Decant supernatant, click tubes, vortex, re-suspend in 10 ml wash MEMA
24. Centrifuge tubes for 10 min at 2000 rpm, 4°C
25. Decant supernatant, click tubes, vortex, re-suspend in 10 ml wash MEMA
26. Centrifuge tubes for 10 min at 2000 rpm, 4°C
27. Decant supernatant, click tubes, vortex, re-suspend in 2 ml wash MEMA, pass five times through 3 cc syringe with 21 gauge needle
28. Determine cell concentration : transfere 100  $\mu$ l of cell suspension in vial with 20 ml Isotone II in Coulter vial
29. Vortex tube, transfer 0.5 ml into tube X.5, vortex tube X.5, (10x dilution - 10(5) cells)  
     transfer 0.5 ml into dilution tube X.4, vortex tube X.4 (100x dilution) - 10(4) cells)  
     transfer 0.5 ml to tube X.3, vortex tube X.3 (1000x dilution) - 10(3) cells)  
     transfer 0.5 ml to tube X.2 and vortex. (10000x dilution) - 10(2) cells)  
 Keep tubes on ice.
30. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate).  
Only X.2 should be seeded for control T-tubes.
31. Transfer 200  $\mu$ l of cell suspension (in triplicate) to 20 ml scintillation vial containing 6 ml

cocktail (Eco-Lume)

32. Incubate dishes for 1 week

33. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol.

Stain colonies with 0.05% crystal violet

34. Count colonies. Should be between 25 and 250 colonies per dish to be a valid data point.

↓

~  
↓

	<u># of cells</u> P60s x 3	<u># of cells</u> P100	For MN
For survival → 1 -	x 2	- For HPRT expression	1 -
2 -	x 2		2 -
3 -	x 2		3 -
4 -	x 2		4 -
5 -	x 2		5 -
6 -	x 2		6 -
7 -	x 2		7 -
8 -	x 3, x 2		8 -
9 -	x 3, x 2		9 -
10 -	x 3, x 2		10 -

}  $\approx 0.5 \times 10^5$   
for 24h with G  
& 48h with G,  
 $\approx 3 \mu\text{g/ml}$  Gpt.

Counts of cells recovered from cluster  
(clusters were pre 3 day at 10.5°C)

**TABLE-3**

Expt. #: V79/Cs-137 acute

Date/Time: 4-6-2001

Plate  
form  
per  
P  
(ml)

Tube #	Coulter count for 100 µl cell suspension	Avg. count - 13 (Background)	Cells/ml [Avg. count x 400 ]		
1	5609 5468 5562	5533	2 213 333		0.05
2	6538 6596 6634	6576	2 630 533		0.06
3	6586 6628 6538	6571	2 628 400	5 256 800	0.07
4	6639 6918 6735	6751	2 700 400		0.08
5	6245, 6282, 6161	6216	2 486 533	4 973 067	0.08
6	6956 6709 6923	6863	2 739 867		0.07
7	6367 6431 6532	6430	2 572 133		0.08
8 →	4597 4525 4148	4410	1 764 133	5 292 400	0.11
9	6452 6373 6309	6365	2 546 000		0.08
10	6655 6810 6805	6744	2 697 467		0.07

Bkgr. - 13  
Plate - 500 µl

IN 3 ml ←

Sample #	Count	Total (in 2ml) $\times 10^6$	$\frac{5 \times 10^5 / P 100}{(\text{volume})}$
1	360	7.2	0.069 0.139
2	339	6.7	0.150
3	80	1.6	0.690
4	330	6.6	0.152
5	99	1.9	0.53
6	290	4.8	0.210
7	260	5.2	0.192
8	332	6.6	0.148
9	290	5.8	0.172
10	306	6.1	0.167



Replating  
1st

on  
Apr. 9, 2007

col. 972/4954977

TERMINATING TIME WRITE  
 + 48 min  
 FT  
 125°C

Is that hypoxic effect involved in this experiment?



Experiment: Apr. 2, 2001

Date/Time:

Tube #	Activity Conc. (kBq/ml)	Activity/Cell (mBq/cell)	Survival Uncorrected	Survival Corrected
1	0.000	0.000	1.0000	1.0000
2	0.000	0.000		
3	#DIV/0!	#VALUE!	0.6990	0.6441
4	#DIV/0!	#VALUE!	0.9029	0.8098
5	#DIV/0!	#VALUE!	0.8544	0.8322
6	#DIV/0!	#VALUE!	0.9612	0.8496
7	#DIV/0!	#VALUE!	0.8447	0.7953
8	#DIV/0!	#VALUE!	0.5049	0.6931
9	#DIV/0!	#VALUE!	0.5437	0.5172
10	#DIV/0!	#VALUE!	0.7379	0.6625

	Dose (Gy)	SF
1(contr)	0	1
2(contr)	0	1
3	0.99	0.6990291
4	1.98	0.9029126
5	2.97	0.8543689
6	3.97	0.961165
7	5.95	0.8446602
8	7.93	0.5048544
9	9.91	0.5436893
10	11.99	0.7378641

$$74 \times 10^4 \rightarrow 10 \times (1)$$

$$SF - \quad \begin{matrix} 10 \times (1) \\ \rightarrow 203 \mu\text{l} \\ + \end{matrix}$$

$$\begin{matrix} 69 \\ 79 \end{matrix} \} / 2 = 74 \times 10^4 / \text{ml}$$

(1500 - x) → 203  
dilute 10x  
20.3 μl


$$\begin{matrix} 3 \times 100 = 300 \\ 3 \times 200 = 600 \end{matrix} \} 900$$

$$2000 = P60' = 4 \text{ ml}$$


$$900$$

$$900 \div 12 + 6 = 18$$

$$900 = 18$$
$$\frac{4}{\text{ml}} \Rightarrow 200 / P60'$$



$74 \times 10^4 \rightarrow 10 \times \rightarrow 10 \times \rightarrow$   
take 203 μl + 29.8 ml MEM





30ml S-MEM

① 30ml S-MEM

② 3ml FCS

③ 0.3ml L-glut

④ 0.3ml P/S

roll over  
rolling

10 -

1 ] CONTROLS

2 ]

3 ]

4 ]

5 ]

6 ]

7 ]

8 ]

9 ]

10 ]

Cs-137  
(dose)

10.5°C

11 ] CONTROLS - ROLL THEM

12 ] OVERNIGHT

13

13 - Control - ~~NO DO~~  
AFTER MARQUESTING  
THE CELLS

$$A - 162 \times 10^4 \times 20 \text{ ml} - 32.4$$

$$B - 236 \times 10^4 \times 20 \text{ ml} - 47.2$$

$$C - 168 \times 10^4 \times 20 \text{ ml} - 33.6$$

$$\begin{array}{r} 32.4 \\ 33.6 \\ \hline 66.0 \times 10^6 \end{array} \begin{array}{l} \text{ml} \\ \text{ml} \\ \text{ml} \end{array}$$

60 → 9.09

60 × 10<sup>6</sup> - 15 ml

5.91

STEP

$$\begin{array}{r} 13 \times 4 \times 10^6 \\ 15 \times 4 = \end{array}$$

60 - 15 ml.

$$\begin{array}{r} 206 \\ 322 \\ \hline 528/2 = 3.14 \text{ ml} \end{array}$$

$$\begin{array}{r} 3.14 - 1 \text{ ml} \\ 4.00 - x = 1.27 \end{array}$$

$$1.27 \text{ (} 10^6 \text{ cells) } \rightarrow 10.73 \text{ sheet}$$

↓

ROLL

$$\begin{array}{r} 18.00 \\ 2.10 \\ \hline 20.10 \\ 3.44 \\ \hline 23.54 \\ 3.00 \\ \hline 26.54 \end{array}$$

V79 / C5-137

Propagate  
28-29-III → UNDER HIT  
30-III - split → NO HIT  
31-III - split → NO HIT

~~30-III~~ }  
4-01-2001 -  $30.5 \times 10^6$  / flask (100% confluency)  $\rightarrow$  plate  $10 \times 10^6$  into A, B, C  
 $30.5 \times 10^6$  / flask (100% confluency)

04-2001 - A -  $162 \times 10^4$  / ml  $\times 20$  ml =  $32.4 \times 10^6$  cells  
B -  $236 \times 10^4$  / ml  $\times 20$  ml =  $47.2 \times 10^6$  cells  
C -  $168 \times 10^4$  / ml  $\times 20$  ml =  $33.6 \times 10^6$  cells.

A+C =  $32.4 + 33.6 = 66 \times 10^6$  cells.

$66 \times 10^6$  (pellet)  $\rightarrow$  + 10 ml  $\xrightarrow{\text{mix medium}}$  take 5.09 ml =  $60 \times 10^6$  cells / 9.09 ml  
 $\downarrow$   
+ 5.91 S-HEP  $\rightarrow$   $60 \times 10^6$  / 15 ml  
 $\approx 4 \times 10^6$  / ml

$\frac{3.14 \times 10^6}{\text{ml}} \leftarrow 2/628 \leftarrow \left| \frac{306}{322} \right. \xleftarrow{\text{COUNT}}$

$3.14 \times 10^6 - 4$  ml  
 $4 \times 10^6 \rightarrow x_{\text{ml}} \rightarrow 1.27 \text{ ml} \rightarrow$  into tube + 0.73 ml S-HEP  
Roll tubes in 37°C - 18:00 / ~~4/02/01~~  $\rightarrow$  14/03/01  
DAY


(b) plated  $3 \times 100$  cells / P60 } for PE.  
 $3 \times 200$  cells / P60 }

(c) plated A1-A3 } app.  $5 \times 10^4$  cells for ~~dose~~ exp. with different  
B1-B3 } conc/time cytotoxicity B  
C1-C3 }

## Results of survival

- Ⓐ Control 1 → cells were harvested and plated for SF, s.  
 same cells were used for rolling & the irradiation with Cs-137

Day	Plated	Colonies	
4-2-2001	100	119 108 119	Ⓐ
4-2-2001	200	221 217 237	Ⓑ



- Ⓑ Control for sample 11 → this is control for SF after overnight (hrs) rolling at 37°C

Day	Sample #	Plated cells	Colonies
4-3-2001	11	100	166 100 141

500-204-47	50	50	50	50-204-4
549-99.9	99.9	59.9-59.9	29.9 19.9	9.9
2633	259	176	114	112
2644	288	156	125	90
2591	232	185	116	101
				50

1.030266

1.545406

$1.03 \times 10^6$   
 $8.5 \times 10^6$   
 $0.2 \times 10^6$   
 $0.19 \text{ ml} + \text{HPRT.}$

500	50	50	50	50	50
5.0	5.0	5.0	5.0	5.0	5.0
99.9	99.9	9.9	59.9	89.9	25.9
4793	424	64	312	330	204
4212	420	83	313	350	
4814	367	82	281	340	

↓ V79 - 10  
 ↓ Beh - 10

$7.4 \times 10^9 \text{ /ml}$   
 Mean count.

$\frac{1.932667}{\text{ml}}$        $\frac{1.614667}{\text{ml}}$   
 $\times 2$                        $\times 2$   
 ↓                              ↓  
 3865333                  3229333

$1.1 \times 10^6 \text{ /ml} \rightarrow$   
 $10 \times \text{①}$   
 $0.1 \text{ ml} - 10^5$   
 $\rightarrow 10 \times \text{②} \rightarrow 10^4$   
 $10 \times \text{③} \rightarrow 10^3 \rightarrow 1000 \text{ /ml} \rightarrow 2.2 \text{ ml} = 200$   
 $10 \times \text{④} \rightarrow 10^2 \rightarrow 100 \text{ /ml}$   
 $\rightarrow 3 \times 100$

CoulterSurvival

Experiment: 0  
 Date/Time: Apr. 2, 2001

Tube #	Coulter count			Average Cells/ml		Hemocytometer Count in Grid			
	1st	2nd	3rd			1st	2nd	3rd	4th
						<i>total # of cells/tube (in 2 ml)</i>			
1	5609	5468	5562	5546	2213333	4.42 × 10 <sup>6</sup>			
2	6538	6596	6634	6589	2630533	5.26 × 10 <sup>6</sup>			
3	6586	6628	6538	6584	2628400	5.24 × 10 <sup>6</sup>			
4	6639	6918	6735	6764	2700400	5.40 × 10 <sup>6</sup>			
5	6245	6282	6161	6229	2486533	~ 5.0 × 10 <sup>6</sup>			
6	6956	6709	6923	6863	2739867	5.48 × 10 <sup>6</sup>			
7	6367	6431	6532	6443	2572133	5.14 × 10 <sup>6</sup>			
IN 3 ml ← 8	4597	4525	4148	4423	1764133 × 3ml	5.28 × 10 <sup>6</sup>			
9	6452	6373	6309	6378	2546000	5.08 × 10 <sup>6</sup>			
10	6655	6810	6805	6757	2697467	5.38 × 10 <sup>6</sup>			

**NOTE:**  
 It seems that cells are propagated after 16 hrs

4 × 10<sup>6</sup> seeded → (4.42) 5.0 ÷ 5.4 × 10<sup>6</sup> - cells are propagated  
 H.I. ...

Tube #	Predicted # Cells Seeded	Actual # Cells Seeded	Colony count 1st	Colony count 2nd	Colony count 3rd	Average	PE (%)	SF Uncorrected	SF Corrected
1	200	221	29	42	31	34	14.176	1.00	1.0000
2	200	263	43	32	29				
3	200	263	21	31	20	24	9.131	0.6990	0.6441
4	200	270	27	37	29	31	11.480	0.9029	0.8098
5	200	249	35	25	28	29	11.797	0.8544	0.8322
6	200	274	43	34	22	33	12.044	0.9612	0.8496
7	200	257	17	29	41	29	11.275	0.8447	0.7953
8	200	176	19	23	10	17	9.825	0.5049	0.6931
9	200	255	19	24	13	19	7.332	0.5437	0.5172
10	200	270	31	19	26	25	9.392	0.7379	0.6625

1. Colonies were fixed & stained on Friday, April 13, 2001. e.g. one week from the day the cells were plated.
2. Very low PE for tube 1 & 2 → question - is that possible that V79 cells can die when they are in cluster for 72h at 10.5°C?

TABLE-4

STAET.

Expt #: V79, C5-132  
 acuk, 100%clus. Date: April, 2, 2001

PE

Controls {

Tube.dilution	Colony 1	Colony 2	Colony 3	Avg Colony	SF <i>NON CORRECTED</i>
1.2/10000	29	42	31	} 34	1
2.2/10000	43	32	29		
3.2/10000	21	31	20	24	0.70
4.2/10000	27	37	29	31	0.90
5.2/10000	35	25	28	29	0.85
6.2/10000	45	34	22	33	0.96
7.2/10000	17	29	41	29	0.85
8.2/10000	19	23	10	17	0.51
9.2/10000	19	24	13	19	0.54
10.2/10000	31	19	26	25	0.74

For tub 8,9,10 also 1000 x diluted cells were plated for survival.

					PE
8.3/1000	142	131	140	138	14%
9.3/1000	151	172	150	158	16%
10.3/1000	154	152	164	157	16%

Note:

For all samples (e.g 8.3, 9.3, & 10.3) number of colonies is not very accurate. Some of colonies were attached each other. Therefore, counts are rather conservative.

Sample #	Count $\times 10^4$ in 2.5ml	Total # of cells $10^6$ in 2.5ml	Volume for $5 \times 10^5$ cells p 100's
1	51x5	6.4	0.195
2	60x5	7.5	0.167
3	45x5	9.4	0.133
4	69x5	8.6	0.145
5	73x5	9.1	0.137
6	73x5	9.1	0.137
7	98x5	12.3	0.102
8	52x5	6.5	0.192
9	80x5	10.0	0.125
10	80x5	10.0	0.125

Replating  
Time  
April 12, 2007



57\_155.

16. T.S. Matney, T.V. Nguyen, T.H. Connor, W.J. Dana and J.C. Theiss, Genotoxic classification of anticancer drugs. *Teratogen. Carcinogen. Mutagen. 5* (1985), pp. 319-328: BIOTECHNOBASE
17. F. Rosner and H. Grunwald, Multiple myeloma terminating in acute leukemia. Report of 12 cases and review of the literature. *Am. J. Med. 57* (1974), pp. 927-939.
18. C.A. Coltman, jr. and D.O. Dixon, Second malignancies complicating Hodgkin's disease: a Southwest Oncology Group 10-year followup. *Cancer Treat. Rep. 66* (1982), pp. 1023-1033.
19. D. Schmahl, M. Habs, M. Lorenz and I. Wagner, Occurrence of second tumors in man after anticancer drug treatment. *Cancer Treat. Rev. 9* (1982), pp. 167-194.
20. A. Niepiński, D.W. Blyant, L. Davison, B. Young, J. Heddle, D.R. McCalla, G. Douglass and K. Michalko, Comparison of three assays for genetic effects of antineoplastic drugs on cancer patients and their nurses. *Environ. Mol. Mutagen. 15* (1990), pp. 83-92. BIOTECHNOBASE

cells in tissue culture. The effect of X-irradiation and clone production with HeLa  
S.A. 41 (1987), pp. 432-437.

3 µg/ml - CH10

5 µg/ml - PL61 ⇒ CHO clone + single copy of prismatic cr. neo + gpt

Set	24h	48h	72h	
Dish	A	B	C	
1	← 1 µg/ml →			2 µl/dish
2	← 3 µg/ml →			6 µl/dish
3	← 6 µg/ml →			12 µl/dish

-1 day - plated  $\approx 5 \times 10^4$   
0 day = Apr. 3 P60

STOCK  
crT B = 2 µg/µl

# V79 COLONY FORMING ASSAY

Experiment Name : 100cluster, no labeling);  
Investigator: M.Lenarczyk  
4/2/2001

Exp.  
Date:

1. Set the rocker-roller at 37°C incubator with 5% CO<sub>2</sub>, set the Coulter Counter, wash cells (from two 150 cm<sup>2</sup> flask, sub-cultured 1:2, 24h before) with PBS, trypsinize cells, each resuspend in 9 ml MEMB, pool, pass 5x through 3 cc syringe with 21 gauge needle, count the cells by transferring 100 ul in Coulter cup containing 20 ml isotone (Coulter balanced electrolyte solution)

2. Dilute to ~4,000,000 cells/ml in MEMB [Actual count :  $4 \times 10^6$  cells/ml]  
Flask A+C = (32.4 + 33.6) x 10<sup>6</sup> cells = 66 x 10<sup>6</sup> cells  
66 x 10<sup>6</sup> cells + 10ml S-MEM → take 9.09 ml (~ 60 x 10<sup>6</sup>) + 5.91 ml S-MEM = 60 x 10<sup>6</sup> / 15ml  
After checking →  $\frac{306}{322}$  → 628/2 →  $\frac{314 \times 10^4}{ml}$ , → take 1.27 ml + 0.73 ml S-MEM (~ 4 x 10<sup>6</sup>/ml)

3. Transfer 1 ml of cell suspension into ten 12 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall

4. Keep the tubes in the roller overnight at 37°C, 5% CO<sub>2</sub> **Date/Time: April 02, 2001 / 18:00**

5. After overnight incubation period (65 hrs) → should be 3 hrs + 12 hrs = 15 hrs. remove tubes and centrifuge at 2000 rpm at 4°C for 10 min (precooled centrifuge). **Date/Time: April 3, 2001 : 10:30**

~~6. Remove buckets from centrifuge and carefully remove 150 µl of supernatant and place in pre-labeled gamma-tube.~~

Make a table with # of samples.

✓  
✓  
2  
3

7. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

8. Centrifuge tubes for 10 min at 2000 rpm, 4°C

9. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

10. Centrifuge tubes for 10 min at 2000 rpm, 4°C

11. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA

12. Centrifuge tubes for 10 min at 2000 rpm, 4°C

13. Decant supernatant, click tubes, vortex, resuspend in 7 ml of MEMA

14. Centrifuge tubes for 10 min at 2000 rpm, 4°C

15. Decant supernatant, click tubes, vortex, transfer the cell suspension in polypropylene microcentrifuge tubes with attached caps (Helena Plastics, 400 ul) using 200 ul pipet tips


16. Again add 200 ul ice cold MEMA, resuspend and transfer the cell suspensions in the same polypropylene microcentrifuge tubes (Total volume ~400 ul)

17. Centrifuge tubes for 5 min at 1000 rpm, 4°C

18. Transfer tubes at 10.5°C for 72 h. **Date/Time: April 3, 2001 /**

→ sample 11 - control for DE + HPRT background + MN(?) →  
AFTER ROLL OVERNIGHT -  
# of cells after roll - is →

100% cluste



10 9.3

10.2 (= 10.8)

157, ~158

$$= 3 \cdot \frac{473}{3} = 3$$

$$\bar{x} = \frac{158}{255} \times 100$$

154

152

164

$$\frac{470}{3} \cdot 3 = 156$$

$$\bar{x} = \frac{157}{270} \times 100 =$$

~~284~~

2546 - 100

158-

158

~~158~~ 158

$$\frac{255-100}{158-x}$$