


## V79 HPRT Mutants/Cell Hypoxic versus Aerobic Clusters



9/28/99

## V79 COLONY FORMING ASSAY

Experiment Name :  $^{137}\text{Cs}$  toxicity (acute, cluster, suspension);

Exp. #: 1;

Experiment performed by: A. Bishayee

Date: 09/06/99

1. Set the rocker-roller at 37°C incubator with 5% CO<sub>2</sub>, set the Coulter Counter, wash cells (from two 150 cm<sup>2</sup> flask, subcultured 1:2, 24h before) with PBS, trypsinize cells, resuspend in 7 ml MEMB for each flask, pool, vortex, pass five times through 3 cc syringe with 21 gauge needle, perform cell count by transferring 100 ul in Coulter cup containing 20 ml isotone (Coulter balanced electrolyte solution)
2. Dilute to ~4,000,000 cells/ml in MEMB (final volume 11 ml) | Actual count: 3,997,333 cells/ $\mu\text{l}$  cells/ml
3. Transfer 1 ml of cell suspension into ten 14 ml tubes (Falcon plastic test tube, 17x100 mm) labeled 1-10 both on cap and wall
4. Roll the tubes for 16 h at 37°C, 5% CO<sub>2</sub> Date/Time: 09/06/99; 4-00 P.M.
5. After ~16 h incubation period, remove tubes, add 8 ml wash MEMA, vortex and centrifuge at 2000 rpm at 4°C for 10 min (precooled centrifuge). Date/Time: 09/07/99; 10-00 A.M.
6. Decant supernatant, click tubes, vortex, resuspend in 3 ml wash MEMA
7. Centrifuge tubes for 10 min at 2000 rpm, 4°C
8. Decant supernatant, click tubes, resuspend in 200 ul ice cold MEMA, transfer the cell suspension in polypropylene microcentrifuge tubes with attached caps (Helena Plastics, 400 ul) using pipet tips
9. Again add 200 ul ice cold MEMA, resuspend and transfer the cell suspensions in the same polypropylene microcentrifuge tubes (Total volume ~400 ul)
10. Centrifuge tubes for 5 min at 1000 rpm, 4°C
11. Transfer tubes at 10°C for 72 h. Date/Time: 09/07/99; 12-00 noon
12. After 72 h, for tubes 1-5, carefully remove the supernatant, resuspend the pellet in 400  $\mu\text{l}$  MEMA and place all tubes on the perforated plate of Rainin pipet tip box containing ice (to maintain ~10.5°C)

13. The tubes were irradiated using Mark I irradiator ( $^{137}\text{Cs}$  gamma-ray), two tube (one tube for pellet and one for the suspension) at a time for a single dose-point, while placing onto a Rainin pipet tip box containing ice as per the Table below

Tube #	Total Dose (R)	Dose rate (Rad/min)	Time (min)	Attenuat.
1	0	0	0	0
2	0	0	0	0
3	300	97.3	3.08	X-10
4	600	739.8	0.81	X-0
5	1200	739.8	1.62	X-0
6	0	0	0	0
7	0	0	0	0
8	300	97.3	3.08	X-10
9	600	739.8	0.81	X-0
10	1200	739.8	1.62	X-0

Resusp

left in pellet

14. After irradiation, carefully remove the supernatant from the top for tubes 6-10, resuspend pellet in 200 ul wash MEMA and transfer the content from all tubes to ten 14 ml tubes (Falcon plastic test tube, 17x100 mm, labeled 1-10 both on cap and wall) containing 10 ml wash MEMA by using pasteur pipet
15. Again add 200 ul wash MEMA in microcentrifuge tubes, resuspend and transfer the cell suspensions in 14 ml tubes
16. Centrifuge the tubes for 10 min at 2000 rpm, 4°C (precooled centrifuge)
17. Labeling and preparation of dilution tubes and colony dishes
  - load 60 mm petri dishes with 4 ml MEMA,
  - load T-tubes with 4.5 ml MEMA and label them 1.2, 1.3, 1.4, 1.5; 2.2, 2.3, 2.4, 2.5; X.2, X.3, X.4, X.5 etc.
18. Decant supernatant, click tubes, vortex, resuspend in 10 ml wash MEMA
19. Centrifuge tubes for 10 min at 2000 rpm, 4°C
20. Decant supernatant, click tubes, vortex, resuspend in 2 ml wash MEMA, pass five times through 3 cc syringe with 21 gauge needle.

21. Determine cell concentration by transferring 100  $\mu$ l to Coulter cup
  22. Vortex tube, transfer 0.5 ml into dilution tube X.5, vortex tube X.5 and transfer 0.5 ml to tube X.4, vortex tube X.4 and transfer 0.5 ml to tube X.3 and vortex tube X.3 and transfer 0.5 ml to tube X.2. Keep tubes on ice.
  23. Transfer 1 ml from dilution tubes into dishes labeled X.2, X.3, X.4 (in triplicate). Only X.2 should be seeded for control T-tubes.
- 
24. Incubate petridishes for 1 week
  25. After 1 week, wash colonies 3 times with normal (1X) saline, and 2 times with methanol. Stain colonies with 0.05% crystal violet.
  26. Count colonies. There must be between 25 and 250 colonies for the flask to be a valid data point.

Date	Day (Gy)	hemo	Coalt (u)	Av	cells/ml	$10^6$	Rel incr.	X4 Total/1000	#doublings	Total do
9/10/90										
	0 1	70x5 = 350	8282 5987	7701	3.08x10 <sup>6</sup> /ml	0.32	1.00	1.2x10 <sup>7</sup>	3.6	
	0 2		7355 7404	7380	2.95	0.34	1.00	1.2	3.6	
	3 3		7962 7852	7907	3.16	0.32	1.08	1.3	3.6	
	6 4		5530 5520	5525	2.21	0.45	0.73	0.88	3.1	
	12 5		5341 5494	5513	2.21	0.45	0.73	0.88	3.1	
	No Resp 0 6		8409 8165	8287	3.31	0.30	1.00	1.3	3.6	
	0 7		8760 8755	8758	3.50	0.29	1.00	1.4	3.8	
	3 8		7748 5395	7139	2.86	0.35	0.81	1.1	3.5	
	6 7		5235 5987	6004	2.40	0.42	0.71	0.96	3.3	
	12 10	62x5 = 310	6424 6240	5850	2.34	0.43	0.70	0.94	3.2	
	BK		5539 2233							
9/15			#s 1 & 9 are contaminated. (S)							
9/16			#s 5 & 6 "							
9/17			switched to Gibco MEM + 5% FBS (Gen 1567A)							
all counted d7	2	5327 5110	4678 3497	6607	4644 ± 1334	1.86x10 <sup>6</sup> /ml	0.54	X4 7.4x10 <sup>6</sup>	2.9	6.5
	3	8254 4838	8354 4756	4875 3112	5698 ± 2125	2.28	0.44	9.12	3.2	6.8
	4	7797 4053	6087 4060	5855 7204	5843 ± 1556	2.34	0.43	9.36	3.2	6.3
	7	8752 5577	4007 5504	4013 5640	4748 ± 909	1.90	0.53	7.6	2.9	6.7
	8	8597 4435	4575 6666	4683 6113	5740 ± 1633	2.34	0.43	9.36	3.2	6.7
	10	9784 4144	6314 6415	4753 6735	6358 ± 1968	2.54	0.39	10.16	3.3	6.5
	BK		11 3 11 6							
9/20		Harvest, dilute $\frac{1}{10}$ ( $\frac{0.1}{1.0ml}$ ) count by hemocytometer								
d10										
	2	221/4 235/4 (adj)	5.7x10 <sup>6</sup> /ml	2.3x10 <sup>7</sup>	cells/1000	4.5		11		
	3	139/4 144/4	3.5x10 <sup>6</sup>	1.4x10 <sup>7</sup>	"	3.8		10.6		
	4	121/4 109/4	2.9x10 <sup>6</sup>	1.2x10 <sup>7</sup>	"	3.6		9.9		
	7	179/4 177/4	4.5x10 <sup>6</sup>	1.8x10 <sup>7</sup>	"	4.2		10.9		
	8	175/4 196/4	4.6x10 <sup>6</sup>	1.8x10 <sup>7</sup>	"	4.2		10.9		
	10	175/4 181/4	4.5x10 <sup>6</sup>	1.8x10 <sup>7</sup>	"	4.2		10.9		

9/20/99  
d10

2	5.7	$\times 10^6$	$\xrightarrow{0.035}$	$2 \times 10^5$	$\xrightarrow{1.000 \times 5}$	
3	3.5		$\xrightarrow{0.057}$			
4	2.9		$\xrightarrow{0.069}$			
7	4.5		$\xrightarrow{0.044}$	0.69		
8	4.6		$\xrightarrow{0.043}$			
10	4.5		$\xrightarrow{0.044}$			
2	5.7	$\times 10^6$	$\xrightarrow[1]{0.0082}$	$5 \times 10^4$	$\xrightarrow[7]{0.0046}$	33
3	3.5	$\times 10^6$	$\xrightarrow[1]{0.014}$		$\xrightarrow[7]{0.0046}$	33
4	2.9	$\times 10^6$	$\xrightarrow[1]{0.017}$		$\rightarrow$	
7	4.5	$\times 10^6$	$\xrightarrow[1]{0.011}$		$\rightarrow$	
8	4.6	$\times 10^6$	$\xrightarrow[1]{0.0100}$		$\rightarrow$	
10	4.5	$\times 10^6$	$\xrightarrow[1]{0.011}$		$\rightarrow$	

3 ml  
66 | P35 x 3

GTG 2.5mg FW 167.2 Add 25ml MEMS to vial = 100  $\mu$ g/ml  
 $100 \mu\text{g/ml} \xrightarrow[325 \text{ ml MEMS}]{5 \text{ ml}}$  167  $\mu\text{g/ml} = 10 \mu\text{M}$

Plate counts

9/27/99 = day 17

Sample	PBS's			Av	std	F100s					Av	std	percent pl
2	80	74	70	74.7	5.03	3	4	7	10	6	6.0	2.74	2.65E-5
3	65	53	64	60.7	6.66	17	18	20	16	20	18.2	1.79	9.89E-5
4	64	56	65	61.7	4.93	28	35	25	27	27	28.4	3.85	15.2E-5
7	64	67	53	61.3	7.37	11	16	24	15	11	15.4	5.34	8.29E-5
8	74	86	69	76.3	8.74	5	9	8	7	11	8.0	2.24	3.46E-5
10	65	65	65	65.0	0	10	11	10	14	17	12.4	3.05	6.30E-5
		overall	av	66.61	8.40								