

TABLE-3

Expt. #: 3

Date/Time: 10/09/00

| Tube # | Coulter count for 100 ul cell suspension | Avg. count | Cells/ml [Avg. count x 4000] | pCi/cell [uCi/ml x 10 ⁶ Cells/ml] | kBq/Cluster [pCi/cell x 148] |
|--------|--|------------|------------------------------------|--|---------------------------------|
| 1 | 577, 592, 563 | | | | |
| 2 | 611, 607, 653 | | | | |
| 3 | 581, 593, 617 | 597 | 2388000 | 0.012 | |
| 4 | 633, 645, 619 | 632 | 2529333 | 0.018 | |
| 5 | 511, 537, 549 | 532 | 2129333 | 0.048 0.079 | |
| 6 | 544, 562, 573 | 559 | 2238666 | 0.094 | |
| 7 | 666, 672, 693 | 677 | 2708000 | 0.204 | |
| 8 | 601, 572, 633 | | | | |
| 9 | 511, 529, 541 | | | | |
| 10 | 532, 555, 562 | 549 | 2198666 | 0.010 | |
| 11 | 513, 549, 562 | 541 | 2165333 | 0.024 | |
| 12 | 562, 539, 547 | 549 | 2197333 | 0.042 | |
| 13 | 560, 542, 522 | 541 | 2165333 | 0.104 | |
| 14 | 680, 669, 671 | 673 | 2693333 | 0.195 | |