

TABLE-3

Expt. # : 2

Date/Time : 04/27/98; 11-00 a.m.

Tube #	Coulter count for 100 ul cell suspension <i>Manometer Select = 50 μl</i>	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]
1	468, 463, 448	459.6	1838666.6	0
2	414, 413, 427	418	1672000	0
3	480, 469, 476	475	1900000	0.00038
4	448, 475, 458	460.3	1841333.3	0.000624
5	432, 424, 447	434.3	1737333.3	0.001640
6	429, 409, 408	415.3	1661333.3	0.003196
7	437, 425, 438	433.3	1733333.3	0.004269
8	375, 366, 394	378.3	1513333.3	0.009310
9	423, 456, 418	432.3	1729333.3	0.01198
10	405, 400, 412	405.6	1622666.6	0.01584