

TABLE-3

Expt. # :

Date/Time :

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]	Km/cell (pa/cell x 148)
1	575, 566, 569				
2	582, 572, 568				
3	595, 582, 601	592	2370666	0.0138	2.05
4	555, 565, 549	556	2225333	0.0663	9.82
5	579, 599, 582	586	2346666	0.1905	28.19
6	611, 623, 607	613	2454666	0.5239	77.54
7	635, 619, 622	625	2501333	0.7827	115.85
8	601, 609, 611				
9	551, 561, 571				
10	569, 557, 547	557	2230666	0.0141	2.08
11	529, 535, 542	535	2141333	0.0665	9.85
12	565, 570, 568	567	2268000	0.1895	28.05
13	582, 592, 579	584	2337333	0.463	68.55
14	575, 562, 555	564	2256000	0.7571	112.05