

TABLE-3

Expt. # : (

Date/Time : 01/01/99; 12-00 noon

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]	mBq/ labeled cell [pCi/cell x 74]
1	518, 550, 569	545	2182666	-	
2	482, 500, 501	494	1977333	-	
3	458, 485, 510	484	1937333	0.00136	0.10064
4	541, 483, 488	504	2016000	0.00634	0.4692
5	531, 525, 523	526	2105333	0.01196	
6	609, 657, 652	639	2557333	0.02158	
7	506, 524, 485	505	2020000	0.06792	
8	432, 440, 428	433	1733333	0.10148	
9	497, 453, 507	485	1942666	0.13687	
10	453, 475, 508	478	1914666	0.27043	