

TABLE-3

Expt. # :

Date/Time :

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 4000]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]
1*	540, 569, 562	557	222800	0.0077
2*	555, 549, 562	555	222133	0.018
3*	515, 529, 507 , 535	526	210533	0.036
4*	475, 489, 467	477	190800	0.046
5*	495, 512, 488	498	199333	0.043
6**	501, 509, 507	505	202266	0.042
7				
8				
9				
10				

mBa/cell
[pCi/cell
x
57]

0.285
 0.666
 1.353
 1.71
 1.57
 1.55

* Processed through serum
 ** washed with HBSS