

TABLE-3

Expt. # : 6

Date/Time : 03/09/98 ; 2-15 p.m.

Tube #	Coulter count for 100 ul cell suspension	Avg. count	Cells/ml [Avg. count x 400]	pCi/cell [uCi/ml x 10 ⁶ Cells/ml]
1	630, 631, 650	637	254800	0
2	777, 791, 796	788	315200	0
3	792, 783, 791	788.6	315466.6	0.1030
4	680, 656, 654	663.3	265333.3	0.2332
5	752, 741, 734	742.3	296933.3	0.3219
6	783, 765, 743	763.6	305466.6	0
7	803, 818, 827	816	326400	0
8	833, 828, 832	831	332400	0.1034
9	750, 784, 764	766	306400	0.2679
10	854, 828, 823	835	334000	0.3254

0.1 = 20 umil
1 = 20
0.1